

VITASEPT[®] WIPES med

Alcohol-based disinfecting wipes

Product properties

- aldehyde- and chlorine-free
- odorant- and colorant-free
- quick exposure time (30 sec)
- doesn't leave residue
- soft fabric-like material

Vitasept[®] Wipes med is an effective broad-spectrum disinfectant for alcohol-resistant medical devices (procedure tables, surgery trolleys, trays, handpieces in dental care, etc.). Soft wipes are permeated with disinfectant and do not leave residues on surfaces after use.

Ingredients

Isopropanol 60%, ethanol 15%

Spectrum of efficacy

Bacteria: *E.coli*, *Salmonella spp.*, *Staphylococcus spp*, *MRSA*, *Streptococcus spp*, *Listeria spp*, *Enterococcus faecalis*, *Corynebacterium ammoniagenes*, *Mycoplasma spp*, *Pseudomonas aeruginosa*, *Enterobacter aerogenes*, *Mycobacterium*, *Clostridium spp* etc

Viruses: HBV, HCV, HIV, Adeno, Rota, Influenza virus etc

Fungi: *Aspergillus niger*, *Candida albicans* etc

Exposure time

Surface disinfection 30sec

Usage

Pull the wipes out until perforated spot according to the amount necessary and if needed support with other hand in order to tear away. All surfaces should be wiped with Vitasept[®] Wipes med, make sure that the surface is covered with disinfectant. Throw the used wipe away. In order to avoid the drying of wipes, close the flip.

Warning

Highly flammable liquid and vapour. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed.

Physical-chemical properties

Physical appearance: *damp cotton wipes*

Color: *colourless*

Scent: *distinctive (alcohol)*

pH: *~7*

Density: *0,85g/cm³ (20°C)*

Storage

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight.

Utilization

Unused product residues and empty packages should be given to licensed waste management.

Conformity to standards

EN 13624, EN 13727, EN 14348, prEN 14476

Other test methods

RKI/DVV

Package

120 wipes in a can, dimensions of a wipe 20 x 20cm
450 wipes in a bucket, dimensions of a wipe 18 x 20cm

