

# GLOBALACID<sup>®</sup>SF

**Amine-based cleaning and disinfection agent for surfaces**

*concentrate*

## Product properties

- cost-efficient, wide spectrum of efficacy
- aldehyde-free
- great cleaning properties
- materialfriendly (metal, plastic, plexi)
- effective also against Tbc, HIV, HBV, Noro virus
- neutral scent
- suitable for use in cleaning machines

## Ingredients

quaternary ammonium compounds 6%, amine 5,5%, polyglycol ether, guanidine, complexing agent.

## Spectrum of efficacy

Unique, patented amine-based disinfectant is virucidal, bactericidal (including *Mycobacterium terrae*) and fungicidal.

## Usage

Globalacid<sup>®</sup>SF is used for the cleaning and disinfecting of various water-resistant surfaces. For example: doors, windowsills, cupboards, floors, bed frames, kitchen inventory, switches etc.

Globalacid<sup>®</sup>SF is available as a concentrate. Make the working solutions according to the exposure time needed, look at the table.

In order to make a working solution out of Globalacid<sup>®</sup>SF, use water that is cool or at room temperature. Once the exposure time has passed, all surfaces that come in contact with food must be rinsed.

**NB! Avoid contact with aldehyde-based products!** Accordingly, surfaces previously treated with aldehyde-based products must be cleaned of aldehyde residues.

**Working solution concentrations and exposure times**

<b>General disinfection, preliminary cleaned environment</b>	0,25% (2,5ml/l)	15min
<b>Infection cases, dirty environment</b>	1,25% (12,5ml/l)	15min
<b>Sporicidal efficacy</b>	3% (30ml/l)	60min

## Warning

Harmful if swallowed. Causes severe skin burns and eye damage. Very toxic to aquatic life. Wear protective gloves/ protective clothing/ eye protection/ face protection.

## Physical-chemical properties

Physical appearance: *colourless, clear liquid*

Scent: *distinctive to amine*

pH: *10,4 (20 °C)*

Solubility: *completely soluble in water*

Density: *0,99 g/cm<sup>3</sup> (20 °C)*

## Storage

See the best before date on the package. Store at room temperature, keep away from direct sunlight and avoid freezing to assure the quality of the product.

## Compliance to standards

EN 1276, EN 14476, EN 13727, EN 13624, EN 14348, EN 14561, EN 14562

## Other test methods

DVG, DGHM/ÖGHMP, DVV/RKI

## Package

1L bottle

5L canister