

GLOBALACID®SF med

Amine-based cleaning and disinfection agent

concentrate

Product properties

- cost-efficient, wide spectrum of efficacy
- aldehyde-free
- great cleaning properties
- materialfriendly (metal, plastic, plexi)
- effective against Tbc, HIV, HBV, Noro virus
- neutral scent

Ingredients

quaternary ammonium compounds 6%, amine 5,5%, polyglycol ether, guanidine, complexing agent.

Spectrum of efficacy

Product is virucidal (including HBV, HCV, HIV, Herpes Simplex, Vaccinia, Noro and Papova virus etc), bactericidal (including *Mycobacterium terrae*, MRSA, *P.aeruginosa*, *E.coli* etc), sporicidal (*B.subtilis*) and fungicidal.

Usage

Globacid®SF med is used for the cleaning and disinfecting of various water-resistant medical devices (patient chairs, stretchers, bed frames etc) and surfaces (doors, windowsills, cupboards, floors, bed frames, kitchen inventory, switches etc). Suitable for cleaning and disinfecting operating room inventory.

Globacid®SF med is available as a concentrate. Make the working solution of clean water and concentrate according to the exposure time needed (look at the table). In order to make a working solution out of Globacid®SF med, use water that is cool or at room temperature. Once the exposure time has passed, all surfaces that come in contact with food must be rinsed.

NB! Avoid contact with aldehyde-based products. Accordingly, surfaces previously treated with aldehyde-based products must be cleaned from aldehyde residues.

Working solution concentration and exposure times

General disinfection, preliminary cleaned environment	0,25% (2,5ml/l)	15min
Infection cases, dirty environment	1,25% (12,5ml/l)	15min
Sporicidal efficacy	3% (30ml/l)	60min

Warning

Harmful if swallowed. Causes severe skin burns and eye damage. Very toxic to aquatic life. Wear protective gloves/ protective clothing/ eye protection/ face protection.

Physical-chemical properties

Physical appearance: *colourless, clear liquid*
 Scent: *distinctive to amine*
 pH: *10,4 (20 °C)*
 Solubility: *completely soluble in water*
 Density: *0,99 g/cm³ (20 °C)*

Storage

See the best before date on the package. Store at temperatures 5-30°C, keep away from direct sunlight.

Utilization

The working solution is ecologically safe, all active ingredients are biodegradable in flush water. For recycling, rinse package with water.

Compliance to standards

EN 1276, EN 14476, EN 13727, EN 13624, EN 13704, EN 14348, EN 14561, EN 14562

Other test methods

DVG, DGHM/ÖGHMP, RKI/DVV

Also conforms to requirements in food industry, beauty services etc.

Package

1L bottle
5L canister

