

GLOBACID[®] ID

ready-to-use solution

Amine-based cleaning and disinfection agent for instruments

ready-to-use solution

Product properties

- wide spectrum of efficacy
- highly effective properties for removing organic debris
- aldehyde, phenol and chlorine free product
- for the pre- and final disinfection of instruments and accessories made from metal and artificial materials, as well as for the disinfection of all endoscopes
- suitable for use in instrument and ultrasound baths as well as endoscope washing machines

Ingredients

quaternary ammonium compounds 18%, amine 14%, nonionic surfactants, complexing agent, corrosion inhibitors.

Spectrum of efficacy and exposure times

Product has a virucidal (HBV, HCV, HIV, papova virus etc), bactericidal (MRSA, *E.coli*, *E.hirae*, *P.aeruginosa*, *Mycobacterium* etc) and fungicidal (*C.albicans* etc) efficacy.

Usage

Spray the product onto disassembled instruments or submerge them into a lid-equipped tub filled with the working solution. Make sure all the instruments are covered with the solution. Once the exposure time has passed, wash or rinse the instruments under running water (at least tap water quality) or with sterile water, as needed and dry with lint-free towel or disposable tissues (instruments can't be left out to dry in the air). It's not recommended to keep the instruments/accessories in the solution longer than the exposure time.

Note! Do not immerse dental handpieces and turbines in the solution. Process them with a rapid disinfectant (like Vitasept[®] Wipes med) by wiping their outer surfaces.

Exposure time

5 min

NB! Avoid contact with aldehyde-based products!

Accordingly, surfaces must be cleaned of aldehyde residues.

When working with chemicals, use protective gloves, eye and/or face protection. Avoid release to the environment.

Physical-chemical properties

Physical appearance: *light-yellow liquid*

Scent: *distinctive to amine*

pH: *ca 10,9*

Solubility: *completely soluble in water*

Density: *1,011 g/cm³ (20°C)*

Storage

See the best before date on the package. Store at temperatures 5-30°C, keep away from direct sunlight.

Utilization

The ready-to-use solution is ecologically safe, all active ingredients are biodegradable in flush water. For recycling, rinse package with water.

Compliance to standards

EN 1276, EN 1650, EN 13727, EN 14348, EN 14562, EN 13624, EN 14561, EN 14476

Other test methods

DVG, DVV/RKI, DGHM/ÖGHMP, Karl Storz endoscope

Package

500ml spray bottle

5L canister



Goodpoint Chemicals OÜ

<http://www.goodpointchemicals.com>

Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia

Contact +372 6626 511 | Fax: +372 6626 522

info@goodpointchemicals.com